

## Hoowaki MicroGrip™


Anti-migration technology for implanted medical devices.



Less worry, better care

Our Hoowaki MicroGrip™ technology prevents device migration and lets you focus on what's important: giving the best care.





200x enlargement of Hoowaki Micropgrip Film Surface

## **Medical Applications**

GI and pulmonary stents

Barrier films

Cardiac repair devices

Cardiac assist devices

Balloons

Hernia meshes


Laparoscopic and robotic grasper and stapler jaws

## Simplifies Approval

- Formed during existing device manufacturing
- · No added chemicals, coatings, or materials
- · Made of same material as device
- Never sticky, abrasive, or aggressive
- Offers a simple route for device qualification; often 510(k)

## **Proven Performance**

Esophageal Stent Pull Against Synthetic Tissue Wet With Saline



Anti-Migration: Up to **5x** resistance to movement



Sold By:

PCH15006